



Logic™ Cart

ACCESSORY FOR PURIFIER®
LOGIC™ BIOSAFETY CABINETS

Logic

The Logic Cart was specifically designed to provide convenient storage of supplies commonly used in biological safety cabinets. Compartments hold and dispense pipettes, gloves, lab tissues and other small supplies. Its compact size allows it to fit easily under Purifier Logic Biosafety Cabinets mounted on Labconco Base Stands. Casters and handles allow the cart to be easily moved to within reach while the user works at the biological safety cabinet.



- Catalog #8022000.
- Three right side compartments, 3.5" w x 17.0" d x 4.5" h (9 x 43 x 51 cm), accommodate pipettes and other supplies.
- Top left side compartment with front panel cutouts accommodates boxes of wipes and gloves. The hinged panel may be pivoted down for loading.
- Drawer, 14.5" w x 15.5" d x 4" h (37 x 39 x 10 cm), holds a variety of miscellaneous small instruments and supplies.
- Open storage area in bottom, 15.8" w x 16.0" d x 12.0" h (40 x 41 x 30 cm), holds larger equipment and supplies.
- Work surface and bottom are constructed of corrosion-resistant white high density polyethylene that is easy to clean. Remaining components are constructed of corrosion-resistant white ABS.
- Compact dimensions allow easy fit under Purifier Logic Biosafety Cabinets mounted on Labconco Base Stands. Work surface is 25.0" w x 19.0" d (64 x 48 cm). Cabinet is 21.5" w x 17.4" d (55 x 44 cm). Overall height is 27.0" (68.6 cm).
- Work surface has handle cutouts in front and right hand side for easy access.
- 2" diameter (5 cm) casters provide mobility.
- Shipping weight 50 lbs. (23 kg)



Logic Carts are shown stored beneath two 6' Purifier Logic Biosafety Cabinets on Hydraulic Lift Base Stands.



Protecting your
laboratory environment

LABCONCO®

Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
800-821-5525 or 816-333-8811
FAX 816-363-0130
www.labconco.com
labconco@labconco.com



© 2007 by Labconco Corporation.
Printed in the U.S.A. Product design
subject to change without notice.